

Product datasheet for **MR203089L3V**

Nudt18 (BC036718) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nudt18 (BC036718) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nudt18
Synonyms:	MGC38179
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC036718
ORF Size:	738 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203089).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC036718 , AAH36718
RefSeq Size:	1630 bp
RefSeq ORF:	740 bp
Locus ID:	213484
Cytogenetics:	14 D2



[View online »](#)

Gene Summary:

Mediates the hydrolysis of oxidized nucleoside diphosphate derivatives. Hydrolyzes 8-oxo-7,8-dihydroguanine (8-oxo-Gua)-containing deoxyribo- and ribonucleoside diphosphates to the monophosphates. Hydrolyzes 8-oxo-dGDP and 8-oxo-GDP with the same efficiencies. Hydrolyzes also 8-OH-dADP and 2-OH-dADP. Exhibited no or minimal hydrolysis activity against 8-oxo-dGTP, 8-oxo-GTP, dGTP, GTP, dGDP and GDP. Probably removes oxidized guanine nucleotides from both the DNA and RNA precursor pools (By similarity).
[UniProtKB/Swiss-Prot Function]