

Product datasheet for **MR203105L2V**

Stc1 (NM_009285) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Stc1 (NM_009285) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Stc1
Synonyms:	Stc
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_009285
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203105).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009285.3 , NP_033311.3
RefSeq Size:	3743 bp
RefSeq ORF:	744 bp
Locus ID:	20855
UniProt ID:	O55183
Cytogenetics:	14 D2
Gene Summary:	Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia. [UniProtKB/Swiss-Prot Function]



[View online »](#)