

Product datasheet for **MR203106L4V**

Gzmb (NM_013542) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gzmb (NM_013542) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gzmb
Synonyms:	AI553453; CCP; CCP-1/C1; CCP-1/C11; CCP1; Ctl; Ctla; Ctla-1; Ctla1; GZB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013542
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203106).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013542.2
RefSeq Size:	1418 bp
RefSeq ORF:	744 bp
Locus ID:	14939
UniProt ID:	P04187
Cytogenetics:	14 28.19 cM



[View online »](#)

Gene Summary:

This gene encodes a member of the granzyme subfamily of proteins, part of the peptidase S1 family of serine proteases. The encoded preproprotein is secreted by natural killer (NK) cells and cytotoxic T lymphocytes (CTLs) and proteolytically processed to generate the active protease, which induces target cell apoptosis. This protein also processes cytokines and degrades extracellular matrix proteins, and these roles are implicated in chronic inflammation and wound healing. Mice lacking a functional copy of this gene exhibit impaired immune cell-mediated cytotoxicity. [provided by RefSeq, Sep 2016]