

Product datasheet for MR203109

Rasl11b (NM_026878) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rasl11b (NM_026878) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rasl11b
Synonyms:	1190017B18Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203109 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGCCTCATCCAGAACATGTGTACCATCGCAGAGTACCCCGCCCAGGCAGCACCGCTGCCGACTGCT
GCCTGGGGGCTGCGGGCCGCGGCTCGTTAAGATCGCCGTGGTAGGGCCAGCGGCTGGCAAGACCGC
TTTGGTGGTCCGATTCTCACTAAACGATTCATCGGGGACTACGAACGAAATGCAGGTAATCTCTACACA
AGACAAGTTCACATAGAAGGAGAAACGCTGGCTATTCAGGTCCAGGACACTCCGGGCATCCAGGTGCACG
AGAATGGCTTGAGTTGCAGCGAGCAGTTGAATCGCTGCATTCGCTGGGCAGATGCCGTCTCCTCGTGTT
CTCCATCAGGACCACAAGAGCTACGAATCATTAGCCAGCTCCACCAGCACGTGCAGCAGCTGCACCCG
GGCACTCGGCTGCCTGTGGTGTGGTGGCTAACAAAGCTGACCTCCTACACGTAAGCAGGTGGACCCCTC
AGCTTGGCCTGCAACTGGCCAGCATGCTCGGCTGCTCCTTCTATGAAGTATCTGTCAAGTAAATACAA
CGACGTCTACAATGCTTTCCATGTCTGTGCAAGAAGTGAGCCCAAGCAGCAGCCAGCAGCACCCCC
GAGAAACGCAGAACTTCCTCATCCCCAGACCCAAAGTCGCCCAACATGCAAGACCTGAAGAGGAGTTTA
AGCAGGCTCTGTCTGCCAAAGTGAGGACTGTACAGTCCGGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203109 protein sequence
 Red=Cloning site Green=Tags(s)

MRLIQNMCTIAEYPAPGSTAADCCGLGAAGRRLVKIAVVGASGVGKTALVVRFLTKRFIGDYERNAGNLYT
 RQVHIEGETLAIQVQDTPGIQVHENGLSCSEQLNRCIRWADAVLVFSITDHKSYELISQLHQHVQQLHP
 GTRLPVVLVANKADLLHVKQVDPQLGLQLASMLGCSFYEVSVSENYNDVYNAFHVLCKEVSPKQPSSSTP
 EKRRTSLIPRPKSPNMQDLKRRFKQALS AKVRTVTSV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_026878

ORF Size: 744 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026878.1](#), [NP_081154.1](#)

RefSeq Size: 1860 bp

RefSeq ORF: 744 bp


Locus ID: 68939

UniProt ID: [Q922H7](#)

Cytogenetics: 5 C3.3

MW: 27.4 kDa

Product images:



Circular map for MR203109