

Product datasheet for **MR203111L3V**

Ms4a6d (NM_026835) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Ms4a6d
Synonyms:	1110058E16Rik; AI447446; Ms4a11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026835
ORF Size:	741 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR203111).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_026835.2](#), [NP_081111.1](#)

RefSeq Size: 1401 bp

RefSeq ORF: 744 bp

Locus ID: 68774

UniProt ID: [Q99N07](#)

Cytogenetics: 19 A



Gene Summary:

May be involved in signal transduction as a component of a multimeric receptor complex.
[UniProtKB/Swiss-Prot Function]