

Product datasheet for **MR203136L3V**

Mrps27 (BC049149) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mrps27 (BC049149) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mrps27
Synonyms:	2610028H14Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC049149
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203136).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC049149 , AAH49149
RefSeq Size:	1833 bp
RefSeq ORF:	746 bp
Locus ID:	218506
Cytogenetics:	13 D1



[View online »](#)

Gene Summary:

RNA-binding component of the mitochondrial small ribosomal subunit (mt-SSU) that plays a role in mitochondrial protein synthesis. Stimulates mitochondrial mRNA translation of subunit components of the mitochondrial electron transport chain. Binds to the mitochondrial 12S rRNA (12S mt-rRNA) and tRNA(Glu). Overexpressed in hepatocellular carcinoma tissues compared with adjacent non-tumoral liver tissues.[UniProtKB/Swiss-Prot Function]