

Product datasheet for **MR203149L4V**

2510009E07Rik (NM_001001881) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	2510009E07Rik (NM_001001881) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	2510009E07Rik
Synonyms:	A1194952
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001001881
ORF Size:	747 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203149).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001001881.1
RefSeq Size:	5164 bp
RefSeq ORF:	747 bp
Locus ID:	72190
UniProt ID:	Q6GQU0
Cytogenetics:	16 B1
Gene Summary:	May play a role in neuronal and neurobehavioral development.[UniProtKB/Swiss-Prot Function]



[View online »](#)