

Product datasheet for **MR203195L3V**

C230093N12Rik (BC031157) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: C230093N12Rik

Synonyms: MGC38572

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: BC031157

ORF Size: 750 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR203195).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [BC031157](#), [AAH31157](#)

RefSeq Size: 3503 bp

RefSeq ORF: 752 bp

Locus ID: 98952

Cytogenetics: 2 B

Gene Summary: May play a role in estrogen action.[UniProtKB/Swiss-Prot Function]

