

Product datasheet for **MR203197A1V**

Mouse Ccn5 (NM_016873) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ccn5 (NM_016873) AAV Particle
Tag: Myc-DDK
Symbol: Ccn5
Synonyms: Ccn5; Crgr4; Ctgfl; Rcop1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR203197 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGGGGCAACCCACTGATCCATCTTCTGGCCATTCCTTCCTCTGCATTCTCTCAATGGTGTATGCC
AGCTGTGCCAGCACCCTGTGCCTGTCTTGACACCACCCAGTGCCACCGGGGTACCCCTGGTGCT
GGATGGCTGTGGCTGCTGTCGAGTGTGTGCACGGAGGCTGGGGAGTCTGCGACCACCTGCATGCTGC
AACCCAGCCAGGGCCTGGTTTGTGACGCTGGGGCAGGCCCAAGTGCCCGTGGTGTGTGCCTCTTCG
AAGAGGATGACGGGAGCTGTGAGGTGAACGGCCGAGGTACCTGGATGGGGAGACCTTTAAACCAATTG
CAGGGTTTTGTGCCGCTGTGATGACGGTGGTTTACCTGCCTGCCGCTGTGCAGTGAGGATGTGCCGGCTG
CCCAGCTGGGACTGCCACGCCCCAGGAGAATACAGGTGCCAGGAAGGTGCTGCCCGAGTGGGTGTGTG
ACCAGGCAGTGATGCAGCCGGCAATCCAGCCCTCCTCAGCCCAAGGACACCAACTTTCTGCCCTTGTAC
TCCTGCATCTGCCGATGGCCCTGTCCAACTGGAGCACAGCCTGGGGCCCTGCTCAACCACCTGTGGG
TTGGGCATAGCCACCCGAGTATCCAACCAGAACCGATTCTGCCAACTGGAGATCCAGCGTCGCTGTGTC
TGTCCAGACCCTGCCTGGCATCCAGGAGCCACGGCTCATGGAACAGTGCCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203197 protein sequence Red=Cloning site Green=Tags(s) MRGNPLIHLAISFLCILSMVYAQLCPAPCACPWTPPQCPPGVPLVLDGCGCCRVCCARRLGESCDHLHVC NPSQGLVCQPGAGPSGRGVVCLFEEDDGSCEVNGRRYLDGETFKPNCRVLCRCDDGGFTCLPLCSEDVRL PSWDCPRPRRIQVPGRCCPEWVCDQAVMQPAIQPSSAQGHQLSALVTPASADGPCPNWSTAWGPCSTTCG LGIATRVSNNRRCQLEIQRRLLCLSRPCLASRSHGWSNSAF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_016873
ORF Size:	753 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_016873.1, NP_058569.1</u>
RefSeq Size:	1784 bp
RefSeq ORF:	756 bp
Locus ID:	22403
UniProt ID:	<u>Q9Z0G4</u>
Cytogenetics:	2 H3
MW:	27.1 kDa