

Product datasheet for **MR203197L3V**

Wisp2 (NM_016873) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Wisp2 (NM_016873) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wisp2
Synonyms:	Ccn5; Crgr4; Ctgfl; Rcop1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016873
ORF Size:	756 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203197).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016873.1 , NP_058569.1
RefSeq Size:	1784 bp
RefSeq ORF:	756 bp
Locus ID:	22403
UniProt ID:	Q9Z0G4
Cytogenetics:	2 H3
Gene Summary:	May play an important role in modulating bone turnover. Promotes the adhesion of osteoblast cells and inhibits the binding of fibrinogen to integrin receptors. In addition, inhibits osteocalcin production (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)