

Product datasheet for **MR203203A1V**

Mouse Nrip3 (BC072641) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Nrip3 (BC072641) AAV Particle
Tag:	Myc-DDK
Symbol:	Nrip3
Synonyms:	A330103B05Rik; D7H11orf14; ICRFP703B1614Q5.2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203203 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTCAGCTGTCACTCATCCCGGGATCGGCCGCGCGTGGACAGGACTCCTCACCGAGGGCGGCCGCA
AGGAGACAGACATGCGGGAAGCGCGTCCCTGCGGCAGCAGCGCCGATGAAGCAGGCGGTGCAGTTCAT
CCACAAGGACTCAGCCGACCTGCTGCCCTGGACGGCCTCAAAAAGCTGGGCTCGTCCAAGGACACGCAA
CCCCATAACATTCTGCAGAGGCGCCTCATGGAACCAACCTGTCTAAGCTCCGGAGCACTCGTGTCCCTT
GGGCCTCTAAGACCAACAAATTCATCAGGCTAAGTCTGAAGGGCTGAAGAAGTGTGAAGATGAAGACAT
GATTTTAGTGTCTTGCCAGTGTGCTGGGAAGGACGTGAAAGCTTTGGTTGACACTGGCTGTCAATATAAT
CTCATCTCCTCAGCCTGTGTGGACAGATTGGGACTAAAGGACCATGTCAAATCTCACAAGCATGAAGGGG
AAAAGCTTTCTCTACCACGGCATCTCAAAGTAGTTGGCCAGATTGAGCACCTAATGATCACAGTGGGTTC
TCTCCGTCTGGATTGCCAAGCAGCTGTGGTTGATGATAATGAAAAAGCTTATCACTTGGTCTCCAGACT
CTTAGATCCCTGAAGTGCATCATAAACTTGACAAGCACCAGCTGATTGTGGGAAGACAGACAAAGAAG
AAATCCCTTTTGTGGAGACAGTTCTGTGAATGACGACAACACTTCAGAAGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203203 protein sequence Red=Cloning site Green=Tags(s) MAQLSVIPGSAAAWTGLLTEGGRKETDMREAASLRQQRRMKQAVQFIHKDSADLLPLDGLKKLGSSKDTQ PHNILQRRMETNLSKLRSTRVPWASKTNKFNQAKSEGLKKCEDEDMILVSCQCAGKDVKALVDTGCQYN LISSACVDRGLGLKDHVKSHKHEGEKLSLPRHLKVVGQIEHLMITVGSRLDCQAAVVDDNEKSLSLGLQT LRSLKCIINLDKHRLIVGKTDKEEIPFVETVSVNDNTSEA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC072641
ORF Size:	753 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC072641, AAH72641</u>
RefSeq Size:	3947 bp
RefSeq ORF:	755 bp
Locus ID:	78593
Cytogenetics:	7 E3
MW:	27.8 kDa