

## Product datasheet for **MR203215L3V**

### Spata9 (NM\_029343) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Spata9 (NM_029343) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Spata9
Synonyms:	1700030K01Rik; 4930599C08Rik; A930023H06Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029343
ORF Size:	756 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203215).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_029343.3</a> , <a href="#">NP_083619.1</a>
RefSeq Size:	1053 bp
RefSeq ORF:	759 bp
Locus ID:	75571
UniProt ID:	<a href="#">Q9D9R3</a>
Cytogenetics:	13 C1
Gene Summary:	May play at role in testicular development/spermatogenesis and may be an important factor in male infertility.[UniProtKB/Swiss-Prot Function]



[View online »](#)