

Product datasheet for **MR203230L4V**

Elf5 (NM_010125) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Elf5 (NM_010125) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Elf5
Synonyms:	ESE-2; ESE-5; ESE-5.
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010125
ORF Size:	762 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203230).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010125.1 , NP_034255.1
RefSeq Size:	2216 bp
RefSeq ORF:	762 bp
Locus ID:	13711
UniProt ID:	Q8VDK3
Cytogenetics:	2 E2



[View online »](#)

Gene Summary:

Transcriptionally activator that may play a role in regulating the later stages of keratinocytes terminal differentiation (By similarity). Binds to DNA sequences containing the consensus nucleotide core sequence GGA[AT]. Transcriptionally activates the TK promoter.
[UniProtKB/Swiss-Prot Function]