

Product datasheet for **MR203234**

C1qb (NM_009777) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C1qb (NM_009777) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C1qb
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203234 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGACACAGTGGGGTGAGGTCTGGACACACCTGTTACTGCTGCTTCTAGGTTTTCTCCATGTGTCT
GGGCCAAAGCAGCTGCACCGGGCCCCCTGGCATCCCTGGCATCCCTGGGGTCCCTGGGGTTCCTGGCTC
TGATGGCCAACCAGGCACTCCAGGGATAAAGGGGGAGAAAGGGCTCCCTGGACTGGCTGGAGACCTTGTT
GAGTTTGGAGAGAAAGGGGACCCAGGGATCCCTGGGACTCCAGGCAAAGTTGGCCCTAAGGGTCCCGTCG
GCCCTAAGGGTACTCCAGGCCCTCTGGACCCCGCGGTCCCAAAGGCGATTCTGGGGACTACGGGGCTAC
ACAGAAAAGTCGCCTTCTCTGCCCTGAGGACCATCAACAGCCCTTGCACCGAACCAGGTACCTCGCTTC
GAAAAGGTGATCACCAACGCGAACGAGAATATGAGCCACGCAACGCAAGTTCACCTGCAAGGTGCCTG
GCCTCTACTACTTACCTATCATGCCAGCTCCCGGGCAACCTGTGTGTGAATCTCGTTCGTGGCCGCGA
TCGGGACAGCATGCAGAAAGTAGTCACCTTCTGTGACTATGCCCAGAACACCTTCCAGGTGACCACAGGT
GGGGTAGTCTTGAAGCTAGAGCAAGAGGAGTTGTTACCTGCAGGCCACAGACAAGAATCCCTCTGG
GCATTGAGGGTGCCAACAGCATCTTCACTGGCTTTCTGCTTTCCCTGACATGGATGCG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >MR203234 protein sequence
Red=Cloning site Green=Tags(s)

MKTQWGEVWTHLLLLLLGFLHVSWAQSSCTGPPGIPGIPGVPVPGSDGQPGTPIKGEKGLPGLAGDLG
 EFGKEKGDPIGTPGKVGPKGPVGPKGTPGSPGPRGPKGDSGDYGATQKVAFSALRTINSPLRPNQVIRF
 EKVITNANENYEPRNGKFTCKVPLYYFTYHASSRGNLCVNLVRGRDRDSMQKVVTFCDYAQNTFQVTTG
 GVVLKLEQEEVHLQATDKNSLLGIEGANSIFTGFLLPDMDA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_009777

ORF Size: 762 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009777.3](#)

RefSeq Size: 1062 bp

RefSeq ORF: 762 bp

Locus ID: 12260


UniProt ID: [P14106](#)

Cytogenetics: 4 69.05 cM

MW: 26.7 kDa

Gene Summary: C1q associates with the proenzymes C1r and C1s to yield C1, the first component of the serum complement system. The collagen-like regions of C1q interact with the Ca(2+)-dependent C1r(2)C1s(2) proenzyme complex, and efficient activation of C1 takes place on interaction of the globular heads of C1q with the Fc regions of IgG or IgM antibody present in immune complexes.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR203234