

Product datasheet for MR203253A1V

Mouse Ska1 (NM_001164355) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ska1 (NM_001164355) AAV Particle
Tag:	Myc-DDK
Symbol:	Ska1
Synonyms:	2810433K01Rik; AV117428
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203253 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGGACTCAGAGCTAGAAGACTTATGCTCTTATGTTAATGAGAAGATTGGAAATATTAAGAAAATCTGT
CGATAAGAAATCTTGGCCAGGACCCTGCCTTGAAGACCACCTTAAGTAAATAGGAGATGAGATCATTGC
AGTAAATGAGCTTCTAAATAAATTTGAATTGGAATTCAGTATCAAGAACAACTAACAGTTCACCTCAAG
GAGCTCTGTGAATCTCTTAGAGAAGAGTGTGAAGATGTGGAACATCTCAAAGAACACGTCCCTCCCCATC
TGCCTCAAGTGACGGCGACCCAGAGCTTAGTTCACAAGCCCAGCCGGATCCTAAAGAGTCAGCAAAAGC
TGAGGAGCCTGGACTGCCGAAGAAGCCCCCAGAGAGCAGAGAGTTATCAAGGAAATGCAATTTATAACT
ATGGATGAATTCAGTGACGTTCCCTGCGTACATGAAATCCCGTTTAAACATATTGTCAAATTAATGATATTA
TTAAGGAAATTAACAAGGCAGTAGTTAGCAAATAAAGATTATGCATCAACAAAAGCATCTATGAGCTC
TGTGAAGAGAAATCTCTACCAGAGATTTAATGAAGAAACGAAGGACACGAAAGGTCATCACTTCATC
GTGGAAGCTGACATAAAGGAATTCACGGCTTTGAAAGTTGACAAGAGGTTTTATGTGATCATGCACATTT
TGCGGCATTGCCACAGACTGTCAGAGGTCGAGGGGGAGGACTTACCCGATACGTCATCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203253 protein sequence Red=Cloning site Green=Tags(s) MDSELEDLCSYVNEKIGNIKKILSIRNLGQDPALKTTLSKIGDEIIAVNELLNKFELEIQYQEQTNSLK ELCESLREECEDVEHLKEHVPPHLPQVTATQSLVHKPEPDPKESDKAEEPGLPKKPPREQRVIKEMQFIT MDEFSDVPAYMKSRLTYCQINDIIKEINKAVVSKYKIMHQPKASMSSVKRNLYQRFINEETKDTKGHHFI VEADIKEFTALKVDKRFYVIMHILRHCHRLSEVRGGGLTRYVIT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001164355
ORF Size:	765 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001164355.1</u> , <u>NP_001157827.1</u>
RefSeq Size:	2605 bp
RefSeq ORF:	765 bp
Locus ID:	66468
UniProt ID:	<u>Q9CPV1</u>
Cytogenetics:	18 E2
MW:	29.4 kDa