

Product datasheet for **MR203270A1V**

Mouse Psma3 (NM_011184) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Psma3 (NM_011184) AAV Particle
Tag:	Myc-DDK
Symbol:	Psma3
Synonyms:	Lmpc8
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203270 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGCTCCATCGGCACTGGGTATGACTTGTCAGCCTCTACGTTCTCTCCTGATGGAAGAGTTTTCCAAG
TTGAATATGCCATGAAGGCTGTGGAAAACAGTAGTACAGCTATTGGGATCAGATGTAAAGATGGTGTGT
GTTTGGGGTAGAAAACTAGTCCTTTCTAACTTTATGAAGAAGGCTCCAATAAACGCTCTTTTAAATGTT
GATCGACATGTTGGAATGGCAGTTGCAGGTTTGCTGGCAGATGCTCGTTCATTAGCAGACATAGCGAGAG
AAGAAGCATCCAACCTTAGGTCCAACCTCGGCTATAACATTCCTCTAAAACATCTTGACAGAGAGTGCC
CATGTATGTACACGCCTATACACTCTACAGCGCTGTTAGACCTTTTGGCTGCAGTTTCATGTTGGGGTCT
TACAGTGCGAATGATGGTGCACAGCTCTATATGATCGACCCATCAGGTGTTTCATATGGTTACTGGGGCT
GTGCGATTGGCAAAGCCAGGCAGGCTGCAAAGACAGAGATAGAAAAGCTTCAGATGAAGGAAATGACTTG
CCGTGATGTAGTTAAAGAAGTTGCAAAAATAATTTACATAGTACACGATGAAGTTAAAGATAAAGCTTTT
GAACTGGAGCTCAGCTGGGTTGGTGAATTAAGTAAAGGAAGACATGAAATTGTTCCCAAAGACATAAGAG
AGGAAGCAGAGAAATATGCCAAGGAATCTTTGAAGGAAGAAGATGAATCAGATGATGACAATATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>MR203270 protein sequence Red=Cloning site Green=Tags(s) MSSIGTGYDLSASTFSPDGRVFQVEYAMKAVENTSSTAIGIRCKDGVVFGVEKLVLSKLYEEGSNKRLFNV DRHVGMVAGLLADARSLADIAREEASNFRSNFGYNIPLKHLADRVAMYVHAYTLYSAVRPFGCSEMLGS YSANDGAQLYMIDPSGSYGYWGCAIGKARQAAKTEIEKLQMKEMTCRDVVKEVAKIIYIVHDEVKDCAF ELELSWVGELTKGRHEIVPKDIREEAKEYAKESLKEEDESDDDNM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_011184
ORF Size:	765 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_011184.2</u>
RefSeq Size:	1434 bp
RefSeq ORF:	768 bp
Locus ID:	19167
UniProt ID:	<u>O70435</u>
Cytogenetics:	12 C2
MW:	28.4 kDa