

Product datasheet for **MR203299L3V**

2210016L21Rik (BC031162) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	2210016L21Rik (BC031162) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	2210016L21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC031162
ORF Size:	768 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203299).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC031162 , AAH31162
RefSeq Size:	1350 bp
RefSeq ORF:	770 bp
Locus ID:	72357
Cytogenetics:	5 F
Gene Summary:	Plays a role in the regulation of Wnt signaling pathway during early development. [UniProtKB/Swiss-Prot Function]



[View online »](#)