

Product datasheet for **MR203309L4V**

Olig3 (NM_053008) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Olig3 (NM_053008) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Olig3
Synonyms:	Bhlhb7; bHLHe20; oligo3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053008
ORF Size:	822 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203309).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053008.2
RefSeq Size:	2086 bp
RefSeq ORF:	822 bp
Locus ID:	94222
UniProt ID:	Q6PFG8
Cytogenetics:	10 A3
Gene Summary:	May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.[UniProtKB/Swiss-Prot Function]



[View online »](#)