

Product datasheet for **MR203326L3V**

Mfap4 (NM_029568) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mfap4 (NM_029568) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mfap4
Synonyms:	1110007F23Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029568
ORF Size:	774 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203326).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029568.2
RefSeq Size:	1513 bp
RefSeq ORF:	774 bp
Locus ID:	76293
UniProt ID:	Q9D1H9
Cytogenetics:	11 37.96 cM
Gene Summary:	Could be involved in calcium-dependent cell adhesion or intercellular interactions. May contribute to the elastic fiber assembly and/or maintenance.[UniProtKB/Swiss-Prot Function]



[View online »](#)