

Product datasheet for **MR203328A1V**

Mouse Cish (NM_009895) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Cish (NM_009895) AAV Particle
Tag: Myc-DDK
Symbol: Cish
Synonyms: Cis; CIS-1; CIS1; F17; F23; SOCS
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR203328 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTCCTCTGCGTACAGGGATCTTGCTCTTTGCTGGCTGTGGAGCAAATTGGGCGGCCGCTCTGTGGG
CCCAGTCCCTGGAGCTGCCCAGGCCAGCCATGCAGCCCTTACCCACTGGGGCATTCCCAGAGGAAGTGAC
AGAGGAGACCCCTGTCCAGGCAGAGAATGAACCGAAGGTGCTAGACCCCTGAGGGGGATCTGCTGTGCATA
GCCAAGACGTTCTCTACCTTCGGGAATCTGGGTGGTACTGGGGTTCTATTACAGCCAGCGAGGCCCGGC
AGCACCTACAGAAGATGCCGGAGGGTACATTCTAGTTCGAGACAGCACCCACCCAGCTACCTGTTTAC
ACTGTCAGTCAAAACCAACCGTGGCCCCACCAACGTGCGGATCGAGTACGCCGATTCTAGCTTCCGGCTG
GACTCTAACTGCTTGTCAGACCTCGAATCCTGGCCTTCCAGATGTGGTCAGCCTTGTGCAGCACTATG
TGGCCTCCTGTGCAGCTGACACCCGGAGCGACAGCCCGGATCCTGCTCCACCCAGCCCTGCCTATGTC
TAAGCAAGATGCACCTAGTGACTCGGTGCTGCCTATCCCGTGGCTACTGCAGTGCACCTGAACTGGTG
CAGCCCTTTGTGCGCAGGAGCAGTGCCCGCAGCTTACAACATCTGTGCGGCTAGTCATCAACCGTCTGG
TGGCCGACGTGGACTGCTTACCCTGCCCCGGCGTATGGCCGACTACCTCCGACAGTACCCCTTCCAAC
T

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203328 protein sequence Red=Cloning site Green=Tags(s) MVLCVQGSCPLLAVEQIGRRPLWAQSLELPGPAMQPLPTGAFPEEVTEETPVQAENEPKVLDPEDLLCI AKTF SYLRESGWYWSITASEARQHLQKMPEGTFLVRDSTHPSYLF TLSVKTTTRGPTNVRIEYADSSFRL DSNCLSRPRILAFPDVSLVQHYYVASCADTRSDSPDPAPTPALPMSKQDAPSDSVLPVPVATAVHLKLV QPFVRRSSARSLQHLCLVINRLVADVCLPLPRRMADYLRQYPFQL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_009895
ORF Size:	771 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_009895.3</u>
RefSeq Size:	2161 bp
RefSeq ORF:	774 bp
Locus ID:	12700
UniProt ID:	<u>Q62225</u>
Cytogenetics:	9 57.99 cM
MW:	28.5 kDa