

Product datasheet for **MR203341**

Nrbp2 (NM_144847) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nrbp2 (NM_144847) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nrbp2
Synonyms:	BC011468
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203341 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGTGCACTGGAGATGGCTGTACTCGAGATCCAAGCCAACGGGGATACCAGAGTCACAGAAGAGGCCA
TCGCTCGAGCCAGGCACTCACTGAGTGACCCCAACATGCGGGAATTCATCCTCTCCTGCCTGGCCCGGGA
CCCTGCCCGCCGACCCTCAGCCCAACCTCTTCCACCGAGTGCTTTGAGGTGCACTCGCTGAAG
CTGCTGGCAGCTCACTGCTTCATCCAGCACCAGTACCTCATGCCTGAGAATGTGGTAGAGGAAAAGACCA
AGGCCATGGACCTCCATGCAGTTTTGGCTGAGATGCCGAGCCCCATGGACCCCAATGCAGTGGCGGTA
CTCAGAGGTCTCCTTCTTGGAGCTGGACAAATTCCTAGAGGATGTCAGGAACGGGATCTATCCACTGATG
AACTTTGCGGCTGCTCGGCCCTTGGGGCTTCCCGTGTGTTGGCCCCACCCCAAGGAAGCCCAAAAGG
CCAAAACCTCCAACGCCAGAACCTTTGACTCGGAGACCAGGAAGGTGGTCCAGATGCAGTGCAACCTGGA
AAGAAGCGAGGACAAGGCTCGGTGGCACCTTACTCTGCTCTTGGTGCTTGAGGACCGGCTACATCGGCAG
CTGACCTATGATCTGCTCCCAACGGACAGTGCCAGGACCTCGCTGCTGAACTAGTGCAATATGGCTTCC
TGACGAGGATGACAGGACAAAGCTAGCAGCCTTCTGGAGACCCTTTCTCAAGTACCGAGGGACGCA
AGCG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >MR203341 protein sequence
 Red=Cloning site Green=Tags(s)

MCALAMAVLEIQANGDTRVTEEAIARARHSLSDPNMREFILSCLARDPARRPSAHNLLFHRVLFVHSLK
 LLAHCFIQHQYLMPENVVEEKTAMDHLAVLAEMPQPHGPPMQWRYSEVSFLELDKFLEDVRNGIYPLM
 NFAAARPLGLPRVLAPPPEEAQKAKTPTPEPFDSETRKVVQMOCNLERSEDKARWHLTLLLVLEDRLHRQ
 LTYDLLPTDSAQDLAELVHYGFLHEDDRTKLAAFLTTFLKYRGTQA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_144847

ORF Size: 774 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_144847.1](#), [NP_659096.1](#)

RefSeq Size: 2537 bp

RefSeq ORF: 777 bp

Locus ID: 223649


UniProt ID: [Q91V36](#)

Cytogenetics: 15 D3

MW: 29.6 kDa

Gene Summary: May regulate apoptosis of neural progenitor cells during their differentiation.
[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR203341