

Product datasheet for **MR203341L3V**

Nrbp2 (NM_144847) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nrbp2 (NM_144847) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nrbp2
Synonyms:	BC011468
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144847
ORF Size:	777 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144847.1 , NP_659096.1
RefSeq Size:	2537 bp
RefSeq ORF:	777 bp
Locus ID:	223649
UniProt ID:	Q91V36
Cytogenetics:	15 D3
Gene Summary:	May regulate apoptosis of neural progenitor cells during their differentiation. [UniProtKB/Swiss-Prot Function]


[View online »](#)