

Product datasheet for **MR203343L4V**

Aph1c (NM_026674) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Aph1c
Synonyms:	0610008A10Rik; APH-1b; Aph1b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026674
ORF Size:	774 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR203343).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_026674.3](#), [NP_080950.1](#)

RefSeq Size: 5185 bp

RefSeq ORF: 777 bp

Locus ID: 68318

UniProt ID: [Q9DCZ9](#)

Cytogenetics: 9 C



Gene Summary:

Potential subunit of the gamma-secretase complex, an endoprotease complex that catalyzes the intramembrane cleavage of integral proteins such as Notch receptors and APP (amyloid-beta precursor protein).[UniProtKB/Swiss-Prot Function]