

Product datasheet for **MR203354A1V**

Mouse Dnajc9 (NM_134081) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Dnajc9 (NM_134081) AAV Particle
Tag:	Myc-DDK
Symbol:	Dnajc9
Synonyms:	5330419I01Rik; AU020082
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203354 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGGCTGCTGGAGCTGTGCGAGCAGGTGTTGGCACCGCCGACCTCTACCAAGTCTGGGCGTGCGGC
GCGAGGCTTCGGACGGCGAGGTCCGACGCGCTACCACAAAGTGTCCCTGCAAGTGACCCCGACCGAGT
AGAGGAGGACCAGAAAGAGGACGCCACCGCCGCTCCAGATCTTGGGAAGAGTCTACGCGGTTCTGAGT
GACAAGGAACAGAAAGCAGTGTACGATGAACAGGGGACAGTGGACGAAGACTCTGCTGGCCTCAACCAAG
ACCGGGACTGGGATGCATATTGGAGATTACTCTTTAAAAAGATATCGCTAGAGGATATCCAAGCTTTTGA
AAAGACATACAAAGGCTCCGAAGAAGAGCTAAACGATATTAAGCAGGCCTATCTGGACTTCAAGGGCGAC
ATGGACCAGATCATGGAGTCTGTACTTTGTGTACAGTACACAGATGAACCCAGGATAAGAAACATCATT
AAAAAGCTATTGAATCCAAAGAGATTCCAGCCTACAGCGCCTTCGTCAAAGAGTCTAAACAAAAGATGAA
TGCAAGGAAAAGGAGGGCTCAGGAAGAGGCTAAAGAAGCAGAGTTGAGCAGAAAGGAGCTGGGACTGGAA
GAAGGAGTGGATAACTTGAAAGCACTCATCCAGAGCAGACAAAAGGATCGGCAAAGGAAATGGACAGTT
TTCTGGCTCAAATGGAAGCAAAATACTGCAACCTTCCAAGGAGGGAAAAGAACAGCACTCAAGAAGGA
AAAGAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203354 protein sequence Red=Cloning site Green=Tags(s) MGLLELCQVFGTADLYQVLGVRREASDGEVRRGYHKVSLQVHPDRVEEDQKEDATRRFQILGRVYAVLS DKEQKAVYDEQGTVDSDAGLNQDRDWDAYWRLLFKKISLEDIQAFETYKGSEELNDIKQAYLDFKGD MDQIMESVLCVQYTDEPRIRNIIQKAIESKEIPAYSAFVKESQKMNARKRRAQEEAKEAELSRKELGLE EGVDNLKALIQSRQKDRQKEMDSFLAQMEAKYCKPSKGGKRTALKKEKK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_134081
ORF Size:	777 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_134081.2</u>
RefSeq Size:	1631 bp
RefSeq ORF:	780 bp
Locus ID:	108671
UniProt ID:	<u>Q91WN1</u>
Cytogenetics:	14 A3
MW:	30.1 kDa