

Product datasheet for **MR203354L3V**

Dnajc9 (NM_134081) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dnajc9 (NM_134081) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dnajc9
Synonyms:	5330419I01Rik; AU020082
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134081
ORF Size:	780 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203354).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134081.2
RefSeq Size:	1631 bp
RefSeq ORF:	780 bp
Locus ID:	108671
UniProt ID:	Q91WN1
Cytogenetics:	14 A3
Gene Summary:	May play a role as co-chaperone of the Hsp70 family proteins HSPA1A, HSPA1B and HSPA8. [UniProtKB/Swiss-Prot Function]


[View online »](#)