

Product datasheet for **MR203356A1V**

Mouse Lrrc61 (NM_177736) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Lrrc61 (NM_177736) AAV Particle
Tag:	Myc-DDK
Symbol:	Lrrc61
Synonyms:	BC027309
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203356 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGCCTCCTGGAGAGAAGCCAGGAGAGGCTGAGGCGCTGAGCATCACACCCAGCTGCTGAAGTCCC
ACTCAGGAGAATTTGCCCTGGATTCCATCCTGTTGCTGAAGCTTCGGGGTTTAGGGGTGGTGACCTGGG
TTGTTTGGGGGAGTGCCTGAACCTTGAGTGGCTTGACCTATCAGGCAATGCACTTACCCACCTGGGCCCA
CTGGCATCCCTGCGCCAGCTGGCTGTGCTCAACGTCTCCAATAATCGGCTTACAGGGCTGGAGCCACTAG
CAGCTTGTGAGAACCTGCAGAGCCTCAATGCTGCTGGTAACCTGCTGACCACTCCTGGTCAGCTGCAGTG
TCTGGCTGGGCTGCAGGCCCTGGAGCACCTGAGGCTCCGGGACCCTTTGGCCCGGCTCAGCAACCCGCTC
TGTGCAAATGCTTCCTACTGGGCTGTGGTCCGAGAGCTACTGCCTGGCCTCAAAGTCATAGATGGGGAAC
GAGTGAGTGGGCGGGGAGTGAGCTGTACCAAGTTATGCCGGGACCTGGACAGTTCCTTGCGCTCCGGCTC
CAGCCCAGGGCCAGAGCCATCGAGGCCAGCCATGGGTAGAGCCCGGCTACTGGGAGTCTGGCCCATC
CGGAGCAGCTCCATTCTGGAGGAGGCCGCGGACAGTTCCAGGACACACTGCAGGAATGCCTTGACCTGG
ATCGTCAGGCCAGCGATAGCTTGGCCAGGCCAGCAGGCTCTCAGCCCTGCAGAAACCACTCTTCCTT
TGCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>MR203356 protein sequence Red=Cloning site Green=Tags(s) MEPPGEEKPGEAEALSITPQLLKSHSGEFALDSILLKLRGLGVVDLGCLGECLNLEWLDLSGNALTHLGP LASLRQLAVLNVSNNRLTGLEPLAACENLQSLNAAGNLLTPGQLQCLAGLQALEHLRLRDPLARLSNPL CANASYWAVVRELLPGLKVIDGERVSGRGSELYQLCRDLSSLRSGSSPGPRAIEAQPWVEPGYWESWPI RSSSILEEACRQFQDTLQECLDLDRQASDSLAQAQALSPAETTSFVF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_177736
ORF Size:	777 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_177736.2, NP_808404.1</u>
RefSeq Size:	2932 bp
RefSeq ORF:	780 bp
Locus ID:	243371
UniProt ID:	<u>Q8R2R5</u>
Cytogenetics:	6 B2.3
MW:	28 kDa