

Product datasheet for MR203357

Gtdc1 (BC052069) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gtdc1 (BC052069) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gtdc1
Synonyms:	E330008O22Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203357 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCAGCCTTATTGGAGATGCTGCTGCATCTCAGTCTGTGGCGCCTTGTCCCCAGCCAGGGCAGAGGG
TTTCTGAGAAAGTCACCGGAGAACTGTGAATCAAAGTCTGATGAACACCCTGACCTTGATGCTGAACAAGA
AGCCCTGGACAACCCATCAGTCCATAAATCAGGCTCTCTCCTGTGAGTAAAGAAAACCTGCCTCTTGAT
CCCAGCACCTTCTGTGTGGAGCTGAAGATCCGCAAAGACCGCTGCACATTACATGGCCGCACAGGTGGG
AGCATGACAAAGATCCAGAGACGTTCTTAAAGATACTAATGAGTCTAAAGCAGTTAAACTTGAACTCCA
TGTCTCTGTCTCGGAGAGACTTTTACAGACACACCAGATATATTTTCAGAAGCCAAAAAGCACTGGGC
TCTTCAGTCTTGCCTGGGCTACTTACCCCGTAAAGAAGACTATTTCCGCGTACTGTGCATGGCTGACG
TTGTCTCTCAACAGCCAAGCATGAATCTTTGGAGTAGCAATGTTGGAAGCTGTGTACTGTGGCTGCTA
CCCCTCTGTCCCAAAGCCTTGGTCTATCCTGAGATCTTCCAGCTGAGTATCTGTATTCCACACCTGAA
CAGCTTCAAACGGCTCAAGAGTTTCTGCAAGAGGCCAGATATTATAAGAAAACATCTTTATAAGGGTG
AAGTGGCTCCATTTCTTGGCAGCCCTACATGGTAAATTCAGGTCTCTGCTCACAACAGAACCCAGGA
AGATTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203357 protein sequence
Red=Cloning site Green=Tags(s)

MLSLIGDAAASQSVAPCPQPGQRVSEKSPENCESKSDEHPDLDAEQEALDNPSVHKSGSLPVSKENLPLD
 PSTLLCGAEDPQRPLHITWPHRWEHDKDPETFLKILMSLKQLNLFHVSVLGETFTDTPDIFSEAKKALG
 SSVLHWGYLPRKEDYFRVLCMADVIVISTAKHEFFGVAMLEAVYCGCYLCPKALVYPEIFPAEYL YSTPE
 QLSKRLKSFCKRPDIIRKHLKYGEVAPFSWAALHGKFRSLLTTEPREDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: BC052069

ORF Size: 777 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC052069](#), [AAH52069](#)

RefSeq Size: 2146 bp


RefSeq ORF: 779 bp

Locus ID: 227835

Cytogenetics: 2 B

MW: 29.1 kDa

Product images:



Circular map for MR203357