

Product datasheet for MR203428

Leng1 (NM_027203) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Leng1 (NM_027203) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Leng1
Synonyms:	1500001K17Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203428 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATATCTTGCCGAAGAAGAGCTGGCACGTCCGGAACAAGGACAATGTTGCCCGCGTGCCGCGGACG
AGGCCAGGCCCGGGAGGAGGAGAAGGAGCGGGAGCGGAGGGTGCTGCTCGCCAGCAAGAGGCCCGTAC
GGAATTCCTACGGAAGAAAGCCAGGCAACGAACTCAGTGCCTGAGCTGGAAGCAGCAGACCCAGGAGCT
CCAAGTCTGGCCCTGTGGACCTGTTTAGGGAGCTGCTGGAAGAAAGAAAGGAGTCCAGAGGCAATA
AAGAGCATGAGGAAGAGAAGCGACGGGAAAAGGAGAGACAAGAGAAAGCACTAGGTATCCTGACTTACTT
GGGCCAGAGTGCAGCAGAGGCCAACTCAGCCTCCTTGGTACCAGCTCCCCCAGGGCAAAGGACTGC
TGCCCTCCAGGCCAAGCCAGATGAGAAGATCAAGAACCAGCTGGATCCCTTAAAAGAGATGCAGAAGC
ATCTGGCAAAGAAAAGGCACAGCAGTGAAGCCGTCTTCCAGAGAGGAACGGCCTCAGAAACAACGACC
CCGAGAGCCACCTTCCCTAGAAAAGCTCCGTGCTGAGCGTCTTCAGCGGGAAGCAGCCGAGAGAGCTCGA
GCAGAGGCCCTGCTGGCCCGGTACAAGGCCAGGTATCCAGCAAGGCCAGGTGGAAGCAGAAGAGACAG
ATGAACGGAGGCGACGATACAACCTCCAGTTCATCCCAAGCTGGCCCGGCGACCCCGCCAGCAGAACC
AACTCTGCCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >MR203428 protein sequence
 Red=Cloning site Green=Tags(s)


MNILPKKSWHVRNKDNVARVRRDEAQAREEEKERERRVLLAQQEARTFLRKKARQNSVPELEAADPGA
 PSSGPVDLFRELLEERKGVPRGNKEHEEEKRREKERQEKALGILTYLQGSAAEAQTQPPWYQLPPGQKDC
 CPPGSPDEKIKNRLDPLKEMQKHLAKKRHSSESRPSREERPQKQRPREPPSLEKLAERLQREAAERAR
 AEALLARVQGVQSQGGQVEAEETDERRRRYNSQFNPLARRPRQNPPTPAH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_027203

ORF Size: 786 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_027203.2](#)

RefSeq Size: 2213 bp

RefSeq ORF: 786 bp


Locus ID: 69757

UniProt ID: [Q9DB98](#)

Cytogenetics: 7 A1

MW: 30.5 kDa

Product images:



Circular map for MR203428