

Product datasheet for **MR203439A1V**

Mouse Nrm (NM_134122) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Nrm (NM_134122) AAV Particle
Tag:	Myc-DDK
Symbol:	Nrm
Synonyms:	2610307M02Rik; AI429796
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203439 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCACCGGCGCTGCTGCTGGTTCTGCGGCCCTCGCCTCCTTCATCCTGGCCTTTGGCACCGGAGTGG
AATTCGTGCGTTTCACCTCCCTGAGGCCGCTGCTTGGAGGGATCCCGAGTCTGGCGGTCCCGATGCCCG
CCACGGGTGGTTGGCTGCCCTGCAGGACCGGAGCATTCTTGCCTCCTTGGCTTGGGATCTATGTCTGCTG
CTACTGTTTGTGGTACAGCACAGCCTTATGGCGACTGAGGCAGTGAAGGCGTGGACATCCCGGTACTTTG
GCGTCCTTCAGAGGTCGCTCTATGTGGCGTGCACAGCCCTGGCCTTGCAGTTGGTGAATGCGGTACTGGA
GACTACACCCAGAGGCCAGTGTATGGGAGGCTCGGGCTGAACCATGGGCCACCTGGGTACCTCTCCTC
TGCTTTGTGCTCCACGTTGTTTCCTGGCTCCTCATCTTTAGCATCCTTCTGGTCTTTGACTACGCTGAAC
TGATGGGCCTCAAGCAGGTTTATTACCATGTGCTGGGACTGGGCGAGCCCTGTCCCTGAAGTCTCCTCG
GGCTCTGAGACTCTTTTCCACCTTCGACACCCAGTGTGTGGAGCTGCTGACTGTGCTGTGGGTGGTC
CCCACCTTGGGCACCGATCGGCTCCTCCTGGCTCTTCTCTTTACCCTCTACCTGGGCTTGGCTCATGGGC
TTGACCAGCAAGATCTCCGCTACCTTCGGTCTCAGCTACAAAGAAAACCTCACCTTCTCTCCAGACCTCA
GGATGGGAAGCAGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203439 protein sequence Red=Cloning site Green=Tags(s) MAPALLLVPAALASFILAFGTGVEFVRFTSLRPLLGGIPESGGPDARHGWLAAALQDRSILASLAWDLCLL LLFVVQHSLMATEAVKAWTSRYFGVLQRSLYVACTALALQLVMRYWETTPRGPVLWEARAEPWATWVPLL CFVLHVVSLLIFSIILLVFDYAELMGLKQVYYHVLGLGEPLSLKSPRALRLFSLRHPVCVELLTVLWVV PTLGTDRLLLALLFTLYLGLAHGLDQDLRYLRSQQLQRKLHLLSRPQDGEAE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_134122
ORF Size:	786 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_134122.1</u>
RefSeq Size:	1502 bp
RefSeq ORF:	789 bp
Locus ID:	106582
UniProt ID:	<u>Q8VC65</u>
Cytogenetics:	17 B1
MW:	29.4 kDa