

Product datasheet for **MR203450L4V**

Mrrf (NM_026422) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mrrf (NM_026422) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mrrf
Synonyms:	2400002D02Rik; MRFF; MTRRF; RRF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026422
ORF Size:	789 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203450).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026422.1
RefSeq Size:	2251 bp
RefSeq ORF:	789 bp
Locus ID:	67871
UniProt ID:	Q9D6S7
Cytogenetics:	2 B
Gene Summary:	Responsible for the release of ribosomes from messenger RNA at the termination of protein biosynthesis. May increase the efficiency of translation by recycling ribosomes from one round of translation to another (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)