

Product datasheet for **MR203451A1V**

Mouse Ttc33 (NM_026213) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ttc33 (NM_026213) AAV Particle
Tag:	Myc-DDK
Symbol:	Ttc33
Synonyms:	2410099M07Rik; 2900001O04Rik; AI507072; AW538342; Osrf
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203451 representing NM_026213 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTTCCTTTGGGTGGAAGAGGAGAATTGGCGAGAAGGTATCCAAGGCCACGTCCCAGCAGTTTGAAG
CAGAAGCAGCTGATGAGAAGGATGCAGCTGAGAACGAGGATGGGAAGTGGCTTCAGGCCAGCAAACGGAG
GAAGGAAACCTGCAGGAGGGCTGTAACAGAGGAGCCAACAGCTGAAGGATGAAGGGGCTCAGCTGGCT
GAAAACAAACGATACAAGGAGGCAATTCAGAAGTGGGATGAAGCACTACAGCTAACCCAGGTGATGCCA
CCCTTTATGAGATGAAATCACAGTCTCTGTCTTTCATGAAATGTTTCCAGCAGTCCATGCAGCAGA
GATGGCTGTCAAGCGTAACCCACACTCGTGGGAGGCGTGCAAACTTTAGGTCGAGCTCAGCTTGGATTA
GGAGAGATAGTTCTGGCAATTCGAAGCTTTCAATAGCCCTGCACATATATCCGATGAACCTGAACTCT
GGAAGGAAGACCTTTCTTGGGCCAGAAAGCTCCAAGAACAGCAGAAAGTAGCCAGAGGATTGAAAACAA
AGAAATGCCCCCGGAAGGCCCGACCTCTCTCAGGGTCCATCCCGGACTATGACTTTGAAAGCGATGAG
ATTGTTGCTGTTTGTGAGCCGTGGCAGAAAAACAGAAGTCGGTGTCTGCAAATAAAACAATGGTTATTG
TGTCGGCGTCTGGGACTGTAGAGATCGTAAACGAGAAGGAAGAGGGCTCTCCACCCAGATGGCTCCGT
CTTTATCAAAGCCCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR203451 representing NM_026213 Red=Cloning site Green=Tags(s) MASFGWKRRIGEKVSKATSQQFEAAAEKDAENEDGNWLQASKRRKETLQEGCKQRSQQLKDEGAQLA ENKRYKEAIQKWDEALQLTPGDATLYEMKSQVLLSLHEMFPAVHAAEMAVKRNPHSWEAWQTLGRAQLGL GEIVLAIRSFQIALHIYPMNPWLKEDLSWARKLQEQQKVAQRIENKEMPPEGPDLSPGSIPDYDFESDE IVAVCAAVAQKQSVSANKTMVIVSASGTVEIVNEKEEGSSTPDGVSF IKAR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_026213
ORF Size:	786 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026213.3, NP_080489.1</u>
RefSeq Size:	1832 bp
RefSeq ORF:	789 bp
Locus ID:	67515
UniProt ID:	<u>Q9D6K7</u>
Cytogenetics:	15 A1
MW:	29.4 kDa