

Product datasheet for MR203482L3

Spsb2 (NM_013539) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spsb2 (NM_013539) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Spsb2
Synonyms:	AI461677; C9; Grcc; Grcc9; SS; SSB2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203482).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_013539
ORF Size:	795 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_013539.1
RefSeq Size:	1179 bp
RefSeq ORF:	795 bp
Locus ID:	14794
UniProt ID:	O88838
Cytogenetics:	6 59.17 cM
Gene Summary:	This gene encodes a member of the SSB family of proteins that contain a central SPRY (repeats in splA and ryanodine receptors) domain and a C-terminal SOCS (suppressor of cytokine signaling) box. The encoded protein is an adaptor protein in the E3 ubiquitin ligase complex that ubiquitinates inducible nitric oxide synthase and targets it for proteasomal degradation. Mice lacking the encoded protein exhibit lower blood urea nitrogen levels and mild thrombocytopenia due to reduced platelet production. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015]

Product images:



Circular map for MR203482L3