

## Product datasheet for **MR203484L3V**

### **Fgf5 (NM\_010203) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Fgf5 (NM_010203) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fgf5
Synonyms:	angora; Fgf-5; Fgf3a; go; HBGF-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010203
ORF Size:	795 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203484).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010203.2</a>
RefSeq Size:	2553 bp
RefSeq ORF:	795 bp
Locus ID:	14176
UniProt ID:	<a href="#">P15656</a>
Cytogenetics:	5 47.77 cM



[View online »](#)

**Gene Summary:**

This gene encodes a secreted protein that is a member of a family of heparin-binding growth factors. The encoded protein regulates cell proliferation, particularly the growth of hair follicles. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013]