

Product datasheet for **MR203510A1V**

Mouse Dffa (NM_010044) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Dffa (NM_010044) AAV Particle
Tag:	Myc-DDK
Symbol:	Dffa
Synonyms:	A330085O09Rik; DFF35; Dff45; ICAD; ICAD-L; ICAD-S
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR203510 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGCTGTCGCGGGGAGCCAGCGCCCCAGACCCGGACGATGTCCGGCCTCTCAAACCGTGTCTGCTAC
 GCCGCAACCACAGCCGCGATCAGCACGGCGTGGCAGCCTCCAGTCTCGAGGAGCTGAGGAGCAAAGCCTG
 TGAATCCTGGCCATTGATAAGTCCCTGACGCCGATCACCTGGTCTGGCTGAGGACGGGACCATAGTG
 GATGACGATGACTACTTCCTCTGCCTTCCTTCCAATACGAAGTTTGTGGCGTTGGCCTGCAATGAGAAGT
 GGATTTATAATGATTCCGATGGAGGGACGGCTTGGGTTTCCCAAGAGTCCTTTGAGGCAGATGAGCCGGA
 CAGCAGGGCAGGGGTGAAGTGGAAGAATGTGGCCAGGCAGCTGAAAGAAGATCTGTCCAGCATCATCCTG
 CTGTGAGAAGAGGACCTCCAAGCGCTCATCGACATCCCATGTGCAGAGCTGGCTCAGGAACCTGCCCCAA
 GCTGTGCCACTGTCCAGGGGCTGCAGAGCACACTCCAGCAGGTGCTTGACCAGAGAGAGGAAGCCCGCCA
 GTCCAAGCAGCTCCTGGAACCTTACCTCCAGGCCTTGAGAAAGAGGGCAACATCTTGTCCAACCAGAAA
 GAGTCCAAGCTGCCCTCAGTGAAGAGCTGGATGCAGTTGACACAGGCGTCGGCAGAGAGATGGCTTCGG
 AAGTGCTGCTCAGAAGCCAGATCCTTACCACACTGAAGGAGAAGCCTGCCCCAGAGCTGAGTTTATCTAG
 TCAGGATTTGGAGGTGGCAAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>MR203510 protein sequence Red=Cloning site Green=Tags(s) MELSRGASAPDPDDVRPLKPCLLRNHSRDQHGVAASSLEELRSKACELLAIDKSLTPITLVLAEDGTIV DDDYFLCLPSNTKFVALACNEKWIYNDSDGGTAWVSQESFEADEPDSRAGVKWKNVARQLKEDLSSIIL LSEEDLQALIDIPCAELAQELCQSCATVQGLQSTLQQVLDQREEARQSKQLLELYLQALEKEGNILSNQK ESKAALSEELDAVDTGVGREMASEVLLRSQILTTLKEKPAPELSLSSQDLEVGKN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_010044
ORF Size:	795 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_010044.2</u>
RefSeq Size:	1097 bp
RefSeq ORF:	798 bp
Locus ID:	13347
Cytogenetics:	4 E2
MW:	29.2 kDa