

## Product datasheet for **MR203511L2V**

### Comt (NM\_007744) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Comt (NM_007744) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Comt
Synonyms:	Comt1; D16Wsu103e; D330014B15Rik
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_007744
ORF Size:	798 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203511).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007744.1</a>
RefSeq Size:	2128 bp
RefSeq ORF:	798 bp
Locus ID:	12846
UniProt ID:	<a href="#">O88587</a>
Cytogenetics:	16 11.4 cM



[View online »](#)

**Gene Summary:**

Catalyzes the O-methylation, and thereby the inactivation, of catecholamine neurotransmitters and catechol hormones. Also shortens the biological half-lives of certain neuroactive drugs, like L-DOPA, alpha-methyl DOPA and isoproterenol.[UniProtKB/Swiss-Prot Function]