

Product datasheet for **MR203517L3V**

Tmem70 (NM_027415) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem70 (NM_027415) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem70
Synonyms:	1110020A09Rik; 2210416J16Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027415
ORF Size:	762 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203517).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027415.2 , NP_081691.2
RefSeq Size:	1729 bp
RefSeq ORF:	762 bp
Locus ID:	70397
UniProt ID:	Q921N7
Cytogenetics:	1 A3
Gene Summary:	Involved in biogenesis of mitochondrial ATP synthase.[UniProtKB/Swiss-Prot Function]



[View online »](#)