

Product datasheet for **MR203558L3V**

Tssk3 (NM_080442) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tssk3 (NM_080442) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tssk3
Synonyms:	1700014N07Rik; 4930594I21Rik; Stk22c; Stk22d; Tsk-3; Tssk-3; Tssk3b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080442
ORF Size:	804 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203558).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080442.2 , NP_536690.1
RefSeq Size:	1034 bp
RefSeq ORF:	807 bp
Locus ID:	58864
UniProt ID:	Q9D2E1
Cytogenetics:	4 D2.2
Gene Summary:	May be involved in a signaling pathway during male germ cell development or mature sperm function.[UniProtKB/Swiss-Prot Function]



[View online »](#)