

Product datasheet for **MR203561L4V**

Capns1 (NM_009795) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Capns1 (NM_009795) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Capns1
Synonyms:	Capa-4; Capa4; Capn4; Cdps; Css1; D7Ertd146e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009795
ORF Size:	807 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009795.2
RefSeq Size:	1461 bp
RefSeq ORF:	807 bp
Locus ID:	12336
UniProt ID:	O88456
Cytogenetics:	7 17.31 cM
Gene Summary:	Regulatory subunit of the calcium-regulated non-lysosomal thiol-protease which catalyzes limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction. [UniProtKB/Swiss-Prot Function]



[View online »](#)