

## Product datasheet for **MR203569L4V**

### 4933425K02Rik (BC025002) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                                  |                                   |
|----------------------------------|-----------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles              |
| <b>Symbol:</b>                   | 4933425K02Rik                     |
| <b>Synonyms:</b>                 | 1700094G20Rik, 4933424B12Rik      |
| <b>Mammalian Cell Selection:</b> | Puromycin                         |
| <b>Vector:</b>                   | pLenti-C-mGFP-P2A-Puro (PS100093) |
| <b>Tag:</b>                      | mGFP                              |
| <b>ACCN:</b>                     | BC025002                          |
| <b>ORF Size:</b>                 | 807 bp                            |

**ORF Nucleotide Sequence:** The ORF insert of this clone is exactly the same as(MR203569).

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                      |                                                                   |
|----------------------|-------------------------------------------------------------------|
| <b>RefSeq:</b>       | <u><a href="#">BC025002</a></u> , <u><a href="#">AAH25002</a></u> |
| <b>RefSeq Size:</b>  | 1564 bp                                                           |
| <b>RefSeq ORF:</b>   | 809 bp                                                            |
| <b>Locus ID:</b>     | 74446                                                             |
| <b>Cytogenetics:</b> | 3 G3                                                              |



**Gene Summary:**

Sperm-specific Na(+)/H(+) exchanger involved in intracellular pH regulation of spermatozoa (PubMed:20036903). Involved in sperm motility and fertility (PubMed:20036903, PubMed:27010853, PubMed:19409551).[UniProtKB/Swiss-Prot Function]