

## Product datasheet for **MR203573L4V**

### Hvcn1 (NM\_001042489) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hvcn1 (NM_001042489) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hvcn1
Synonyms:	0610039P13Rik; AI450555; BTS; HV1; mVSOP; Vsop.
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001042489
ORF Size:	807 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203573).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001042489.1</a> , <a href="#">NP_001035954.1</a>
RefSeq Size:	2782 bp
RefSeq ORF:	810 bp
Locus ID:	74096
UniProt ID:	<a href="#">Q3U2S8</a>
Cytogenetics:	5 F



[View online »](#)

**Gene Summary:**

Mediates the voltage-dependent proton permeability of excitable membranes. Forms a proton-selective channel through which protons may pass in accordance with their electrochemical gradient. Proton efflux, accompanied by membrane depolarization, facilitates acute production of reactive oxygen species in phagocytosis.[UniProtKB/Swiss-Prot Function]