

## Product datasheet for **MR203580L4V**

### Fibcd1 (NM\_178887) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fibcd1 (NM_178887) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fibcd1
Synonyms:	AI448887
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178887
ORF Size:	810 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203580).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178887.2</a> , <a href="#">NP_849218.1</a>
RefSeq Size:	4764 bp
RefSeq ORF:	1380 bp
Locus ID:	98970
UniProt ID:	<a href="#">A2AV25</a>
Cytogenetics:	2 B



[View online »](#)

**Gene Summary:**

Acetyl group-binding receptor which shows a high-affinity and calcium-dependent binding to acetylated structures such as chitin, some N-acetylated carbohydrates, and amino acids, but not to their non-acetylated counterparts. Can facilitate the endocytosis of acetylated components (By similarity).[UniProtKB/Swiss-Prot Function]