

Product datasheet for **MR203593**

Use1 (NM_001145780) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Use1
Synonyms:	2010315L10Rik; 5730403H22Rik; AV002165; D12; Ed2; p31; Q-snare
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR203593 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCAGGCGGAAGGAGCCTATCATCGCCCGCTGGCGACGTCGCGGTTGGAGTTGAACCTGGTACGGC
TGCTTTGTCGCTGTGAGTCGATGGCGGCAGAAAAACGAGAACCAGGACGAGTGGCGACTGGAGAAGTATGT
GGGTGCCCTGGAAGACATGCTGCAGGCCCTGAAAGTTCAAGCAAGTAAACCGGCCTCAGAAGTATCAGC
GAGTACTCCCGCAAAGTGGACTTTCTGAAGGGCATGCTGCAGGCCGAGAAGTTGACCTCCTCCTGAAA
AAGCTCTAGCTAACCAGTTCTGGCTCCTGGCCGCGTGCCCAACACAGCCAAGGAGCGGGTGCCTGCCAC
CAAGACTGTACATTTGAGTCTAGGGCACGTTACACCAGTGAGATGCGGAGTGAGCTGTTGGGCATGGAA
CCATCTGGAGAATGTGAGGTGGACATGAGGAAACGAGCTGCGAAGGGGTCCAGACCTGCAGATGAGAGGC
AGTCAGCATCCGAGCTAGACCTCGTCTCCAGAGGCACCAGGGCTCCAGGAAAAGCTGGCTGAGGAGAT
GCTAGGCCTGGCCCGCAGCCTCAAGACCAACACACTGGCTGCCAGAGTGTCATCAAGAAGGACAATCAG
ACCCTGTCACATCACTCAAGATGGCAGACCAGAACCTGGAGAAGCTGAAGCTGGAGTCGGAGAGATTGG
AGCAGCATGCACAGAAGTCAGTGAAGTGGCTGCTGTGGGCCATGCTCATTGTCGTGTCTTCGTGTTTCAT
CAGCATGATCCTGTTTATCCGAATCATGCCAGGCTGAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR203593 protein sequence
 Red=Cloning site Green=Tags(s)


MAQAE GAYHRPLATSRLELNLRLLCRCE SMAAEKREPDEWRLEKYVGALEDMLQALKVQASKPASEVIS
 EYSRKVDFLKGMLQAEKLTSSSEKALANQFLAPGRVPTTAKERVPATKTVHLQSRARYTSEMRSELLGME
 PSGECEVDMRKRAAKGSRPADERQSASELDLVLQRHQGLQEKLAEEMGLARSLKTNTLAAQSVIKKDNQ
 TLSHSLKMADQNLEKLKLESERLEQHAQKSVNWLLWAMLIVVCFVFI SMILFIRIMPRLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001145780

ORF Size:

810 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq Size:	1087 bp
RefSeq ORF:	855 bp
Locus ID:	67023
UniProt ID:	Q9CQ56
Cytogenetics:	8 B3.3
MW:	30.6 kDa
Gene Summary:	SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR203593