

Product datasheet for **MR203597**

LOC329575 (BC092391) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	LOC329575
Synonyms:	OTTMUSG00000016325
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR203597 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAAACCAAAATTTAAAGAATTCATACAGATAAAAAGACTCATCAGTATAATCAGTGTAAGAAAGATT
TCACGTGTGCAAAATTTATTTGCAGGCATGAAAGAAGTTGTTCTGCAGAGCAACCCTCTGAGTTTATTCA
ATGTGGTAAAGCCTTTGCATATGAGAGTGGTAGACAAAGGCATCAAATTAACATACTGGAGAGAAACAC
CATGACTGTAACCAATGTGGTAAGGACTTTGGAACATGGAATGCTCCTGCAATACATAAGCGAACACATA
CAGGAGAGAAACCCTATGACTGTAACAATGTGGTAAAGCCTTTGCAAGAAGCTGTCATCTCCGAATACA
TAATCGAACACATGCAGGAGAGAAACAGTATGAATGTAACCAATGTGGTAAAGCTTTTAAAGAAGGAGT
GACCTCCAATACATAAGCGAACACATACAGGAGAGAAACCCTATGAATGTAACCAATGTGGTAAAGCCT
TTGCAAGAAGTGGTGACCAACAAGAACATAAGCGAACACATACAGGAGAGAAACCCTATGAATGTAACCA
ATGTGGTAAAGCTTTTATAAGAAGGAGAGTCTCCAAATACATAAGCGAACACATACAGGAGAGAAACCA
TATGAATGTAACAATGTGGTAAAGCCTTCGCAAGCAGTAGTGACCTCCAAAAACATAAACGGACACATT
CCGGAGAGAAACCCTATGAATGTAACAATGTGGTAAAGCTTTTGACAAAGCAGTCATCTCCGAATACA
TAAGCAAACACATACAGGAGAGAGACCCTATGAATGTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>MR203597 protein sequence
Red=Cloning site Green=Tags(s)

MKPKFKRIHTDKKTHQYNQCKDFTCANFICRHERSCSAEQPSEFIQCGKAFAYESGRQRHQIKHTGEKH
 HDCNQCGKDFGTWNVLQIHKRTHTGEKPYDCKQCGKAFARSCHLRTHNRTHAGEKQYECNQCGKAFKRRS
 DLQIHKRTHTGEKPYECNQCGKAFARSGDQQEHKRTHTGEKPYECNQCGKAFIRRRVLQIHKRTHTGEKP
 YECKQCGKAFASSDLQHKRTHSGEKPYECKQCGKAFQSSHLRTHKQTHTGERPYECN



TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN:

BC092391

ORF Size:

810 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC092391](#), [AAH92391](#)

RefSeq Size: 3536 bp


RefSeq ORF: 812 bp

Locus ID: 329575

Cytogenetics: 2 H4

MW: 31.5 kDa

Product images:



Circular map for MR203597