

Product datasheet for **MR203597L3V**

LOC329575 (BC092391) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	LOC329575
Synonyms:	OTTMUSG00000016325
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC092391
ORF Size:	810 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR203597).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	BC092391 , AAH92391
RefSeq Size:	3536 bp
RefSeq ORF:	812 bp
Locus ID:	329575
Cytogenetics:	2 H4

