

Product datasheet for **MR203615L4V**

Plpp7 (NM_145521) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Plpp7 (NM_145521) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Plpp7
Synonyms:	D830019K17Rik; NET39; Ppapdc3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145521
ORF Size:	816 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203615).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145521.2
RefSeq Size:	1824 bp
RefSeq ORF:	816 bp
Locus ID:	227721
UniProt ID:	Q91WB2
Cytogenetics:	2 B



[View online »](#)

Gene Summary:

Plays a role as negative regulator of myoblast differentiation, in part through effects on MTOR signaling. Has no detectable enzymatic activity. Knockdown in myoblasts strongly promotes differentiation, whereas overexpression represses myogenesis.[UniProtKB/Swiss-Prot Function]