

Product datasheet for **MR203617**

Ube2j2 (NM_021402) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ube2j2 (NM_021402) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ube2j2
Synonyms:	1200007B18Rik; 2400008G19Rik; 5730472G04Rik; AL022923; NCUBE-2; Ubc6; Ubc6p
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR203617 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCAATAACAGCAATAAGAGAGCTCCAACAACAGCAACACAAAGGTTGAAGCAGGACTACCTTCGAA
TTAAGAAAGACCCAGTGCCTTACATCTGTGCAGAGCCCCCTCCCTTCAAATATTCTTGAATGGCATTATGT
TGTCGAGGCCCTGAGATGACTCCTTATGAAGGTGGCTATTATCATGGAAAATAATTTTCCAGAGAA
TTTCCATTTAAGCCTCCTAGTATTTACATGATAACCCCAATGGAAGATTTAAGTGAACACAAGGCTGT
GTCTCTCCATCACAGATTTCCACCCAGATACATGGAATCCAGCGTGGTCCGTCTCCACCATTCTAACAGG
GCTCCTGAGCTTTATGGTGGAGAAAGGCCCCACCTGGGGAGCATAGAGACGTCGGACTTCACGAAAAAA
CAGCTGGCAGCACAGAGTCTAGTGTAAATTTAAAGATAAAAGTCTTTTGTGAATTTTCTGAAGTTG
TGGAGGAAATTAACAGAAAACAGAAAGCACAAGATGAACTAAGTAACAGGCCCGAGAATCTTCTTTGCC
AGACGTGGTTCCTGATGGGGAGCTGCACCGTGGCCAGCATGGGATCCAGCTCCTCAATGGGCACGCCCC
GCTGCTGGCCCCAACCTCGTGGGCTCCACAGGCCAACAGGCATCACGGACTCCTGGCGGGCGGTTGG
CGAACTTGTTTGTATAGTTGGGTTTGCAGCCTTCGCCTACACGGTCAAGTACGTGCTGAGAAGCATAGC
ACAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR203617 protein sequence
 Red=Cloning site Green=Tags(s)

MSNNSNKRAPTTATQRLKQDYLRIKKDPVYPYICAEPLPSNILEWHYVVRGPEMTPYEGGYH GKLI FPRE
 FPFKPPSIYMITPNGRFKCNTRLCLSIDFHPDTWNPAWSVSTILTGLLSFMVEKGPTLGS IETSDFTKK
 QLAAQSLVFNLKDKVFCELFPEVVEEIKQKQKAQDEL SNRPQNLPLPDVVPDGELHRGQHGIQLLN GHAP
 AAGPNLAGLPQANRHHGLLGGALANLFVIVGFAAFAYTVKYVLRISIAQE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_021402

ORF Size: 780 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021402.5](#), [NM_021402.6](#), [NP_067377.4](#)

RefSeq Size: 3520 bp

RefSeq ORF: 816 bp

Locus ID: 140499


UniProt ID: [Q6P073](#)

Cytogenetics: 4 87.66 cM

MW: 29 kDa

Gene Summary: Catalyzes the covalent attachment of ubiquitin to other proteins. Seems to function in the selective degradation of misfolded membrane proteins from the endoplasmic reticulum (ERAD).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR203617