

Product datasheet for **MR203617L4V**

Ube2j2 (NM_021402) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2j2 (NM_021402) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2j2
Synonyms:	1200007B18Rik; 2400008G19Rik; 5730472G04Rik; AL022923; NCUBE-2; Ubc6; Ubc6p
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021402
ORF Size:	780 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203617).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021402.5 , NP_067377.4
RefSeq Size:	3520 bp
RefSeq ORF:	816 bp
Locus ID:	140499
UniProt ID:	Q6P073
Cytogenetics:	4 87.66 cM
Gene Summary:	Catalyzes the covalent attachment of ubiquitin to other proteins. Seems to function in the selective degradation of misfolded membrane proteins from the endoplasmic reticulum (ERAD).[UniProtKB/Swiss-Prot Function]



[View online »](#)