

Product datasheet for **MR203626L4V**

Rem2 (NM_080726) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rem2 (NM_080726) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rem2
Synonyms:	AW411893
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080726
ORF Size:	819 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203626).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080726.2 , NP_542764.1
RefSeq Size:	1887 bp
RefSeq ORF:	1026 bp
Locus ID:	140743
UniProt ID:	Q8VEL9
Cytogenetics:	14 C2
Gene Summary:	Binds GTP saturably and exhibits a low intrinsic rate of GTP hydrolysis.[UniProtKB/Swiss-Prot Function]



[View online »](#)