

## Product datasheet for **MR203632L2V**

### **Spi1 (NM\_011355) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Spi1 (NM_011355) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Spi1
Synonyms:	Dis-1; Dis1; PU.1; Sfpi-1; Sfpi1; Spi-1; Tcfpu1; Tfpu.1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_011355
ORF Size:	816 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203632).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011355.1</a>
RefSeq Size:	1333 bp
RefSeq ORF:	819 bp
Locus ID:	20375
UniProt ID:	<a href="#">P17433</a>
Cytogenetics:	2 50.44 cM



[View online »](#)

**Gene Summary:**

Binds to the PU-box, a purine-rich DNA sequence (5'-GAGGAA-3') that can act as a lymphoid-specific enhancer. This protein is a transcriptional activator that may be specifically involved in the differentiation or activation of macrophages or B-cells. Also binds RNA and may modulate pre-mRNA splicing.[UniProtKB/Swiss-Prot Function]