

## Product datasheet for **MR203641L4V**

### Pitpnb (BC034676) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pitpnb (BC034676) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pitpnb
Synonyms:	PI-TP-beta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC034676
ORF Size:	816 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR203641).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC034676</a> , <a href="#">AAH34676</a>
RefSeq Size:	2681 bp
RefSeq ORF:	818 bp
Locus ID:	56305
Cytogenetics:	5 F
Gene Summary:	This gene encodes a member of the phosphatidylinositol transfer protein family. The encoded protein catalyzes the transfer of phospholipids (phosphatidylinositol and phosphatidylcholine) between membranes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2014]



[View online »](#)